

# SPECIALIZED GERMAN ORTHOPAEDIC




## INSTRUMENTS AND SETS




**Biomed Healthtech (P) Ltd.**  
Marthanda, 207, 2<sup>nd</sup> Floor, Dr. A. B. Road,  
Worli Naka, Mumbai - 400 018. INDIA.  
Tel : +91 - 22 - 2491 1409, 4333 9200  
Fax : +91 - 22 - 2492 1598  
E-mail : [biomedimp@vsnl.com](mailto:biomedimp@vsnl.com)  
Website: [www.biomedhealthtech.com](http://www.biomedhealthtech.com)

# EXTRACTION SET FOR BROKEN SCREWS



**EMERGENCY INSTRUMENTS TO REMOVE  
DAMAGED AND BROKEN BONE SCREWS**

EMERGENCY INSTRUMENTS TO REMOVE DAMAGED AND  
BROKEN BONE SCREWS (UP FROM DIAMETER 1.5 MM TO 7.0 MM  
SCREWS)



SCREW-EXTRACTION-SET

54-09-25-000

complete



54-09-25-0000-31

TRAY(Empty)



RECOMMENDED FOR  
STERILIZATION CONTAINER FOR

54-09-25-000

complete

CONTAINER BOTTOM

305-100

285 x 280 x 100 mm

CONTAINER LID

300-010

285 x 280 mm

**SET - LISTING FOR SET: 54-09-25-000  
SCREW EXTRACTION-SET COMPLETE**



**SET - LISTING FOR SET: 54-09-25-000 SCREW- EXTRACTION - SET COMPLETE**

REF. - NO.:	PRODUCT - DESCRIPTION:	
54-09-25-0000-31	AL-Tray for Screw-Extraction Set-54-09-25-000	-1-Pcs.
52-09-15-1603	Conical Extraction Screw for 1.5 / 2.0 mm Screws	-1-Pcs.
52-09-27-1603	Conical Extraction Screw for 2.7 / 3.5 / 4.0 mm Screws	-1-Pcs.
52-09-45-1603	Conical Extraction Screw for 4.5 / 5.0 / 6.5 / 7.0 mm Screws	-1-Pcs.
52-09-15-1601	Extraction Bolt for 1.5 mm Screws	-1-Pcs.
52-09-20-1601	Extraction Bolt for 2.0 mm Screws	-1-Pcs.
52-09-27-1601	Extraction Bolt for 2.7 mm Screws	-1-Pcs.
52-09-35-1601	Extraction Bolt for 3.5 / 4.0 mm Screws	-1-Pcs.
52-09-45-1601	Extraction Bolt for 4.5 mm Screws	-1-Pcs.
52-09-65-1601	Extraction Bolt for 5.0 / 6.5 / 7.0 mm Screws	-1-Pcs.
52-09-15-1602	Hollow Reamer for 1.5 mm Screws	-1-Pcs.
52-09-20-1602	Hollow Reamer for 2.0 mm Screws	-1-Pcs.
52-09-27-1602	Hollow Reamer for 2.7 mm Screws	-1-Pcs.
52-09-35-1602	Hollow Reamer for 3.5 / 4.0 mm Screws	-1-Pcs.
52-09-45-1602	Hollow Reamer for 4.5 mm Screws	-1-Pcs.
52-09-65-1602	Hollow Reamer for 5.0 / 6.5 / 7.0 mm Screws	-1-Pcs.
52-09-15-1602-03	Spare Reamer (Reamer Tube) for Hollow Reamer 1.5 mm (used for 1.5 mm Screws)	-1-Pcs.
52-09-20-1602-03	Spare Reamer (Reamer Tube) for Hollow Reamer 2.0 mm (used for 2.0 mm Screws)	-1-Pcs.
52-09-27-1602-03	Spare Reamer (Reamer Tube) for Hollow Reamer 2.7 mm (used for 2.7 mm Screws)	-1-Pcs.
52-09-35-1602-03	Spare Reamer (Reamer Tube) for Hollow Reamer 3.5 / 4.0 mm (used for 3.5 / 4.0 mm Screws)	-1-Pcs.
52-09-45-1602-03	Spare Reamer (Reamer Tube) for Hollow Reamer 4.5 mm (used for 4.5 mm Screws)	-1-Pcs.
52-09-65-1602-03	Spare Reamer (Reamer Tube) for Hollow Reamer 5.0 / 6.5 / 7.0 mm (used for 5.0 / 6.5 / 7.0 mm Screws)	-1-Pcs.
30-46-16-1800	Sharp Hook, Length 155 mm	-1-Pcs.
21-03-00-1600	Mini Lexer Gouge, 4 mm	-1-Pcs.
21-03-00-1601	Mini Lexer Gouge, 6 mm	-1-Pcs.
21-03-00-1603	Mini Lexer Gouge, 10 mm	-1-Pcs.
	Screw Holding Forceps 20.5 cm / 8"	-1-Pcs.
30-06-16-1600	T-Handle, Large with AO-Quickcoupling Length 160 mm	-1-Pcs.
54-09-25-000-01	Template for use Broken Screw-Extraction Instruments (Guide Manual for Use)	-1-Pcs.

## SELECTED INSTRUMENTS

### Conical Extraction Screws



52-09-15-1603	Conical Extraction Screw for 1.5 / 2.0 mm Screws
52-09-27-1603	Conical Extraction Screw for 2.7 / 3.5 / 4.0 mm Screws
52-09-45-1603	Conical Extraction Screw for 4.5 / 5.0 / 6.5 / 7.0 mm Screws

### Hollow Reamers complete



52-09-15-1602	Hollow Reamer for 1.5 mm Screws
52-09-20-1602	Hollow Reamer for 2.0 mm Screws
52-09-27-1602	Hollow Reamer for 2.7 mm Screws
52-09-35-1602	Hollow Reamer for 3.5 / 4.0 mm Screws
52-09-45-1602	Hollow Reamer for 4.5 mm Screws
52-09-65-1602	Hollow Reamer for 5.0 / 6.5 / 7.0 mm Screws

### Extraction Bolts



52-09-15-1601	Extraction Bolt for 1.5 mm Screws
52-09-20-1601	Extraction Bolt for 2.0 mm Screws
52-09-27-1601	Extraction Bolt for 2.7 mm Screws



52-09-35-1601	Extraction Bolt for 3.5 / 4.0 mm Screws
52-09-45-1601	Extraction Bolt for 4.5 mm Screws
52-09-65-1601	Extraction Bolt for 5.0 / 6.5 / 7.0 mm Screws

## SELECTED INSTRUMENTS

### Hollow Reamers Tube



52-09-15-1602-03	Speare Reamer (Reamer Tube) for Hollwow Reamer 1.5 mm (used for 1.5 mm Screws)
52-09-20-1602-03	Speare Reamer (Reamer Tube) for Hollwow Reamer 2.0 mm (used for 2.0 mm Screws)
52-09-27-1602-03	Speare Reamer (Reamer Tube) for Hollwow Reamer 2.7 mm (used for 2.7 mm Screws)
52-09-35-1602-03	Speare Reamer (Reamer Tube) for Hollwow Reamer 3.5 / 4.0 mm (used for 3.5 / 4.0 mm Screws)
52-09-45-1602-03	Speare Reamer (Reamer Tube) for Hollwow Reamer 4.5 mm (used for 4.5 mm Screws)
52-09-65-1602-03	Speare Reamer (Reamer Tube) for Hollwow Reamer 5.0 / 6.5 / 7.0 mm (used for 5.0 / 6.5 / 7.0 mm Screws)

### Hollow Reamers Central Bolt



52-09-15-1602-02	Central Bolt used with Hollow Reamer for 1.5 mm
52-09-20-1602-02	Central Bolt used with Hollow Reamer for 2.0 mm
52-09-27-1602-02	Central Bolt used with Hollow Reamer for 2.7 mm
52-09-35-1602-02	Central Bolt used with Hollow Reamer for 3.5 / 4.0 mm
52-09-45-1602-02	Central Bolt used with Hollow Reamer for 4.5 mm
52-09-65-1602-02	Central Bolt used with Hollow Reamer for 5.0 / 6.5 / 7.0 mm

### Sharp Hook



30-46-16-1800

Length 155 mm

## SELECTED INSTRUMENTS

**T-Handle Large with AO-Quickcoupling**



30-06-16-1600

Length 160 mm

**Screw Holding Forceps**



30-01-21-1603

Length 20,5 cm / 8"

**Mini Lexer Gouge**



21-03-00-1600

Length 4 mm

**Mini Lexer Gouge**



21-03-00-1601

Length 6 mm

**Mini Lexer Gouge**



21-03-00-1603

Length 10 mm